

Dedicated OLISA™ Biochips for Multiple GMO Detection and Identification

E. Bufflier,


emmanuel.bufflier@apibio.com

<http://www.apibio.com>

OLISA™ Biochips for Multiple GMO Detection and Identification

- **Principle of OLISA™ Biochip Technology**
- **Example of a dedicated GMO OLISA™ Biochip**
- **OLISA™ : an integrated system**

Principle of OLISA™



Specific Hybridization

Schematic OLISA™ Protocol




DNA
extraction



PCR
amplification





OLISA™
Reactions (automatable)




1-2 hours

*Hands-on
time : 30 min.*



Imaging

OLISA™ GMO Detection and Identification Biochip				
GMO	2030900001	6,2004	T1	
	19/09/2003		D01	
Sample4 : B63-806				
Sample contains corn BT11				
Parameters	Results	Values	Mode	Positions
Controls +	+	32603	DO.A	P1: 1,17
p35S	+	25674	DO.A	P1: 2
676-678-680	-	0	DO.A	P1: 3
DLL25	-	0	DO.A	P1: 4
ADH	+	22563	DO.A	P1: 5
BT176	-	0	DO.A	P1: 6
MON 802-809	-	0	DO.A	P1: 7
MS3-MS6	-	0	DO.A	P1: 8
676-678-680	-	0	DO.A	P1: 9
BT11	+	24987	DO.A	P1: 10



Principle of OLISA™ technology in DNA analysis and Genotyping

- **Hybridization probe : SNP detection resolution**
- **Post amplification detection (PCR, NASBA, RCA, others), biotin-labeled target.**
- **Enzyme-mediated densitometric detection.**

Advantages of OLISA™

- **Miniaturization and Multidetecction**

1 sample → 1 experiment → up to 16 results

- **High-Throughput compatible 96 wells format**


- 500 tests per day manually (500 tests)

- 3000-7000 tests per day using liquid handling robotic system

- **Customizable & Flexible Manufacturing (Easy Array Update)**

- **Routine biochip application compatible with market prices**

Example of a dedicated GMO OLISA™ chip



Probe	Target
S1	tNOS
S2	ADH
S3	Bt176
S4	T25
S5	MON810
S6	GA21
S7	CBH351
S8	BT11
S9	RR
S10	Lectin
S11	p35S

Limit of sensitivity

● **Smaller concentrations detected**


p35S and tNOS : 20 copies

ADH standard maize : 100 copies


Lectin standard soybean : 20 copies

Standard Maize ADH

100 copies




0.1 copies



Intra-test Reproducibility

- 4 concentrations (0 ; 0.01 ; 0.1 and 1%)
- 2 repetitions




Pool of all transgenic Maize lines

Inter-test Reproducibility

- 4 concentrations (0 ; 0.01 ; 0.1 and 1%)

- 2 extractions



- 2 independent tests
(different days)



Test Reliability

● Results with reference method/Results with OLISA

Samples	S35	Nos	Corn B	T25	Mon810	Bt176	Bt11	CBH351	GA21	Lectine	RR
C1	-/-	-/-	+/+	-/-	-/-	-/-	-/-	-/-	-/-	-/-	-/-
C2	+/+	+/+	+/+	-/-	-/-	-/-	-/+	+/+	-/-	-/-	-/-
C3	+/+	+/+	+/+	+/+	+/+	+/+	-/-	-/-	-/-	-/-	-/-
S1	-/-	-/-	-/-	-/-	-/-	-/-	-/-	-/-	-/-	+/+	-/-
S2	+/+	+/+	-/-	-/-	-/-	-/-	-/-	-/-	-/-	+/+	+/+



Dedicated Software Development : Report

- General report for the whole experiment :

	01	02	03	04	05	06	07	08	09	10	11	12
A	A01	A02	A03	A04	A05	A06	A07	A08	A09	A10	A11	A12
B	B01	B02	B03	B04	B05	B06	B07	B08	B09	B10	B11	B12
C	C01	C02	C03	C04	C05	C06	C07	C08	C09	C10	C11	C12
D	D01	D02	D03	D04	D05	D06	D07	D08	D09	D10	D11	D12
E	E01	E02	E03	E04	E05	E06	E07	E08	E09	E10	E11	E12
F	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
G	G01	G02	G03	G04	G05	G06	G07	G08	G09	G10	G11	G12
H	H01	H02	H03	H04	H05	H06	H07	H08	H09	H10	H11	H12


No chips (no strips)



- **Green : GMO-free sample, Red: Sample contains GMO(s)**

Dedicated Software Development : Report

● Dedicated report for each sample :

OLISA™ GMO Detection and Identification Biochip				
GMO	2030900001	6,2004	T1	
19/09/2003		D01		
Sample4 : B63-806				
Sample contains corn BT11				
Parameters	Results	Values	Mode	Positions
Controls +	+	32603	DO.A	P1: 1,17
p35S	+	25674	DO.A	P1: 2
676-678-680	-	0	DO.A	P1: 3
DLL25	-	0	DO.A	P1: 4
ADH	+	22563	DO.A	P1: 5
BT176	-	0	DO.A	P1: 6
MON 802-809	-	0	DO.A	P1: 7
MS3-MS6	-	0	DO.A	P1: 8
676-678-680	-	0	DO.A	P1: 9
BT11	+	24987	DO.A	P1: 10
NK603	-	0	DO.A	P1: 11
TC1507	-	0	DO.A	P1: 12
DBT418	-	0	DO.A	P1: 13
T14-T25	-	0	DO.A	P1: 14
MON 863	-	0	DO.A	P1: 15
tNOS	+	26246	DO.A	P1: 16



Custom development and production

- **All GMO chips,**
- **Authorized/Non authorized GMO chips,**
- **Dedicated Gene-Maker/Seed Developer
GMO chips,**
- **... your GMO chips**

- **Easy Array Update (New GM crops
launching)**

Current customers in GMO field

- Seed developers ...



- Bioanalytical Services Labs ...



Other Applications of Apibio Chips

ILISA : protein based biochips (multidetetection ELISA principle)

Precipitating
Enzymatic
Substrate



Streptavidine-enzyme
Conjugate

Biotinylated
Detection antibody

Antigen


Capture antibody

Principle



Pattern and Result of a
Food Allergen
ILISA™ Biochip

CORRELATION ASSAY



Comparison with a
Reference Method

Contact us ...

- **Location**

**Zone ASTEC, 15 rue des Martyrs,
F-38054 GRENOBLE Cedex 9, France**

Tel. : +33 (0)4 38 78 62 33

 **contact@apibio.com**

- **www.apibio.com**

- **Facilities : 800 m2**



Meet us at our Booth today ! Thank you