




A new versatile oligonucleotide micro-array (OLISA™: OLigo Sorbent Arrays) has been developed by Apibio to simultaneously monitor up to 1572 different genes. This technology uses the 96-well microtiter plate as a platform in which up to 17 different probes can be attached to the bottom of each well. Therefore, OLISA™ technology permits a high throughput multidetection assay. Using synthetic biotin-labeled targets (50-70 mers), OLISA™ was shown to detect a low number of targets (<math>10^7</math> copies per assay) with a dynamic range of 2.5 logs. Both low and high copy number messenger RNAs were detected from 0.2  $\mu$ g of total RNA per array, as demonstrated on a referenced biological model for cadmium toxicity using yeast cells.


## Principle, Material and Methods



## Results




**Fig. 1:** Imaging performance of hybridized biotin-labeled synthetic 50-mer targets (act: even, seo: odd).



**Fig. 2:** Dynamic range of OLISA™ Colorimetric assay (> 2 logs).

**Sensitivity:** 10<sup>7</sup> copies per assay

**Reproducibility:** Coefficients of variations are lower than 10%



**Fig 3:** Gene expression profiling in yeast cells. Simultaneous detection of highly represented mRNA (*yef* > 20 copies / cell)\* and rare mRNA (*seo* ~ 0.1 copy / cell)<sup>□</sup> from 0.2  $\mu$ g of total RNA. In contrast of the 4-fold down-regulated *yef* \* gene, *ynl* <sup>◊</sup> and *seo* <sup>□</sup> genes are 4 and 5-fold, respectively, up-regulated by the cadmium.

## Conclusion

OLISA™ based on 96-well microtiter plate is a flexible and efficient technology that allows in a short time to simultaneously monitor the gene expression profiling for up to 17 different genes per well. Besides its flexibility, OLISA™ technology based on colorimetric assay has major advantages:

- A convenient format on a laboratory standard format that easily can fit most of existing instruments.
- A high sensitivity and a large dynamic range that allow to analyze simultaneously different expressed genes from a small sample.
- A standardized protocol which can be used in relevant domains such as gene expression profiling, drug screening, toxicology and therapeutic monitoring.
- A proprietary densitometry reading system and analysis software.
- A cost-competitive DNA chip technology compatible with high-throughput requirements.